



Bug Buster Kit for detection and removal of head lice

In NHS Choices (July 2010) and NHS Clinical Knowledge Summaries (March 2010), the bevel-edged Bug Buster comb is recommended for the detection and removal of head lice. The Bug Buster Kit contains three brilliant Bug Busters plus the superlative Nit Buster comb. Whatever treatment is used, every family needs a Bug Buster® combing Kit as well. Do not miss out on routine pharmacy sales! Pip Code: 233-1783

Now is the perfect time to take advantage of our discount prices.

Bug Buster Kit: Suggested retail price, off the shelf: **£8.42** inc. VAT
Single BBK pharmacy price: **£4.68** plus VAT 0.59p = £5.27 each inc. delivery

We offer 5% discount on 15 – 49 Kits and 7.5% discount on 50 – 99 Kits

No. of Kits	Cost per BBK	Cost	VAT	Total Price
5	£ 4.68	£ 23.40	£ 2.95	£ 26.35
15	£ 4.45	£ 66.75	£ 8.40	£ 75.15
50	£ 4.33	£216.50	£ 27.00	£ 243.50

Attractive window poster *Bug Buster Kit available here* supplied with your order

Fax/telephone **01908 261501** or email **bugbusters2k@yahoo.co.uk** to order, giving the following details:

Contact Name:.....Date of Order:/...../.....

No. of Kits required:.....Total Price: £.....

Delivery to:.....Branch No:.....

Address:.....

.....Post Code:.....Tel. No:.....

Invoice to (if different from above):.....

.....Post Code:.....Tel. No:.....

Supplier: Community Hygiene Concern, 22, Darin Court, Crownhill, Milton Keynes, MK8 0AD

Not-for-profit reg. charity no: 801371; VAT is calculated using an apportionment method

For same day posting, 1st class, Monday to Friday, place your order by 12 noon

Payment terms: 28 days from invoice date E&O exc. Prices subject to change without notice

Credit/debit card payment appreciated – please ring 01908 561928

Data Sheet

BUG BUSTER® KIT for reliable detection and treatment of head lice without medication, and comfortable removal of nits

Head Lice Device

The Bug Buster Kit contains 3 Bug Buster combs, 1 Nit Buster comb, 1 wide tooth comb, a protective cape, and full instructions on the detection and treatment of head lice using the Bug Busting® wet combing method and subsequent Nit Busting to remove unsightly nits.

The Bug Buster Kit works in combination with ordinary shampoo and ordinary hair conditioner – no additional medicated product is required.

It is *reusable* by a whole family for louse detection and treatment of infestation.

Clinical particulars – therapeutic indications

To detect and remedy head infestation by systematic elimination of head lice, and, if desired, easy nit removal for cosmetic purposes.

Method of administration

Wash hair with usual shampoo, rinse, apply hair conditioner generously, untangle and straighten hair with wide tooth comb. Pass a Bug Buster comb from roots to tips, section by section, until no more lice are found. Clear the comb of lice between strokes, by wiping both sides on paper towel or rinsing. Typically the time required to cover the whole head ranges from 2 minutes for short, straight hair to 30 minutes for long, frizzy hair. Rinse off conditioner and comb wet hair again with a Bug Buster comb and then the Nit Buster comb. Repeat 3 times at half weekly intervals (a total of 4 sessions in 2 weeks). If re-infestation occurs, recognisable by the presence of full-grown lice at session 2, 3 or 4, continue for a further 3 sessions from their detection.

Also used to check the efficacy of any treatment choice 5, 9 and 12/13 days after application.

Marketing authorisation

Registered medical device, ref: CA 006936, 3 May 2001

Sept 2002 Drug Tariff, Part IXA (available on NHS prescription even if not mentioned by name in some editions of the *BNF* or other formularies, because it is an appliance)

Pip code: 233-1783

Net NHS price = £4.68 from October 2011

Support service

A Help Line service on 01908 561928, provides one to one support. A demonstration DVD/video may be purchased from Community Hygiene Concern.

Manufacturer

Produced on a non-profit making basis by:

Community Hygiene Concern, registered charity no: 801371

Unit 22, Darin Court, Crownhill, Milton Keynes, MK8 0AD

Internet: www.chc.org Email: bugbusters2k@yahoo.co.uk

Sales / Fax: 01908 261501

October 2011